

# CT Image Denoising

Clearer Images. Stronger Diagnoses. Brighter Tomorrows.

Team Name


**UNO REVERSE**

Team Members

**Shubham Vishwakarma**

**Jeet Nakrani**

**Dhiraj Nair**



# Index

✓ Problem Statement

✓ Proposed Solution

✓ Innovation / Uniqueness

✓ Tech Stack


✓ Architecture / Workflow

✓ Execution Plan

✓ Scalability & Real World Impact

✓ Anticipated Challenges

✓ Conclusion



# Problem Statement

## Understanding the Problem

In medical imaging, particularly brain CT scans, image quality is often compromised by two prevalent types of noise:

### 1. Poisson Noise:

- Originates from the statistical nature of photon detection during image acquisition.

### 2. Periodic Noise:

- Caused by interference from electronic and mechanical components in the CT scanner.

## Aim of the Solution


The objective is to develop a **software tool** capable of:

- **Reading and displaying CT brain scans** in DICOM format.
- **Designing an algorithm** that reduces Poisson and periodic noise effectively.
- Enhancing image clarity while preserving essential diagnostic information, leading to **more accurate diagnoses**.

# Proposed Solution



We propose a convolutional autoencoder that learns to remove both Poisson noise and Periodic noise from brain CT images by enhancing image clarity, while preserving the important diagnostic features.



The system will read and process CT scans in DICOM format while preserving all essential imaging metadata. Because of this, the system will be directly compatible with real-world clinical workflows.



The framework may also be extended to MRI and ultrasound as well and potentially fit into real-time hospital systems for AI-assisted diagnostics.

## Key Functional Modules



### Encoder

Extracts low-level features and compresses noisy input.



### Bottleneck

Captures critical patterns while discarding noise.



### Decoder

Reconstructs high-quality, denoised CT images.



### Activation Functions


Uses ReLU for layers and Sigmoid for output.



### Loss Function

Mean Squared Error (MSE) ensures pixel-level restoration accuracy.

# Innovation / Uniqueness



# Tech Stack

## Machine Learning & Deep Learning



## Medical Imaging



## Deployment & Infrastructure




## Evaluation & Monitoring





# Architecture / Workflow



# Execution Plan

1

## Data Collection

Gathering and preparing CT brain scans for processing

2

## Model Development

Designing and training a convolutional autoencoder

3

## Evaluation & Testing

Assessing the model's performance and diagnostic clarity

4

## User Interface

Creating a user-friendly interface for interaction

5


## Deployment

Deploying the model on AWS Cloud for accessibility

6


## Monitoring & Maintenance

Ensuring long-term reliability through regular updates





# Scalability & Real World Impact



## Scalability



**Cross-Modality Extension:** The model can be adapted for MRI, PET, and ultrasound imaging with minimal architectural changes.



**Real-Time Integration:** Can be integrated with live CT workflows to enhance the quality of images immediately during the diagnosis.




**Cloud Deployment:** Scalable hosting on platforms like AWS for hospital-wide or cross-institution deployment.




**Modular Architecture:** Can easily be upgraded with future AI advances or additional noise profiles.



**Improved Diagnostic Accuracy:** Cleaner images give doctors greater confidence in detecting tumors, hemorrhages and abnormalities.



**Clinical Workflow Enhancement:** Seamless DICOM handling and user interface reduces radiologist burden and reduces manual filtering.




**Reduced Misdiagnosis Risk:** Denoised images reduce the false-positive/false-negative rates when it comes to assessing critical brain scans.



**Global Applicability:** Can benefit hospitals that lack resources or even telemedicine platforms as a way of providing high-quality imaging periodically.

# Anticipated Challenges



**Data Limitations:** Difficult to get large numbers of clean-labeled DICOM brain scans. Privacy issues and access restrictions in medical datasets.



**High Computational Requirements:** Model building and real-time inference require GPU-based systems. Might hinder the deployment at hospitals where resources are limited.




**Diagnostic Sensitivity:** Over-denoising may occur resulting in subtle abnormality patterns becoming blurred. Very important to balance denoising while preserving medical-based features.



**Regulatory & Ethical Compliance:** Must comply with HIPAA, GDPR, and medical data handling policies. Legal barriers of introducing AI-based tools into clinical settings.




**Clinical Acceptance:** Resistance from radiologists towards AI-based tools. Will require considerable explainability, validation and trust in order to overcome existing resistance



# Conclusion

- Our model successfully **reduced both Poisson and Periodic noise**, producing **much cleaner CT brain scans**.
- While denoising we **maintained the integrity of the scans** and **all relevant diagnostic features were preserved**
- To validate the effectiveness of our approach, we created a **predictive classification model** that assessed how well the **denoised images supported accurate diagnosis** — comparing performance **before and after denoising**.
- The complete solution was **deployed on AWS EC2 instance** where users process the CT scan data by **single and bulk upload**, ensuring **high performance and scalability** in real-world use.

***Because behind every clear scan  
is a life depending on it!***





**Thank You**